Articulating Paper Forceps – Instructions For Use

Cleaning and Disinfection:

- Thoroughly clean after each use and before sterilization using one of the methods listed below:
 - Ultrasonic Cleaning.
 - Automated Washing Equipment.
 - Manual scrubbing in warm soapy water.
- Corrosive cleaners and detergents should not be used due to risk of chemical reactions with the stainless steel.
- Abrasive bristle (wire) brushes are not recommended as the surface finish of the instrument could be damaged.
- Hinged instruments should be lubricated after cleaning and rinsing and before sterilization.
- Thoroughly dry before sterilization.

Sterilization:

- Sterilize in a sterilization pouch or instrument cassette according to the following parameters:
 - Steam Autoclave cycle of 20 minutes @ 121°C.
 - Unsaturated Chemical Vapor cycle of 20 minutes @ 132°C.
 - Dry Heat Sterilization cycle of 90 minutes @ 160°C.
- Hinged instruments must be sterilized in the 'open' position to ensure effective sterilization of the hinge area.
- Do not sterilize along with non-stainless steel metal items.

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